

Intro to Ultrasound - MSK

September 19, 2020
DMU-COM
Thomas Benzoni, DO, AOBEM, FACEP

Declarations

- COI
 - None
- Agency
 - none
- FDA indication deviation
 - None; these are *“radiation emitting devices”*
- Rules of engagement
 - Adult conversation

Uses, now and future

- “Stethoscope of the future”
 - Maybe...
- Answer a specific question
 - All good answers start with a good question.
 - Shock? Which type (tank, pump, pipes)?
 - Dead? Cardiac motion?
 - Pregnant? Where/wellbeing?
 - Access? Seeing is believing...and cannulating.
 - Fluid? Can I put a needle in there?
 - Inflated? Better than X-ray!

Objectives

- Basic physics
- Machine operation
- Probe (transducer) types
- Probe manipulation
- Obtaining an image
- Artifacts
- Reporting/exporting/admin (at session end)

Physics


- Sound wave sent, received (*transducer*)
 - Transducer sends briefly, listens long
- Time received - time sent = distance
- \angle of incidence = \angle of reflection*
 - Source of artifact (anisotropy)
- Reflection occurs at tissue-type interfaces
- Frequency: α data density, $1/\alpha$ penetration*
 - Impacts probe selection

Machine operation


- Knobs (“knobology”)
 - Gain
 - Does NOT change probe; adds “volume” to signal
 - Depth
 - Does NOT change probe
 - Instructs software to ignore signals after certain time interval
 - Doppler
 - Signal approaching = red, receding = blue (convention)
 - M-Mode, etc.

Probe types

- Curved array
 - Curved line source, higher frequency; abd. studies
- Linear array
 - Straight line, highest frequency; vascular, MSK
- Phased array
 - Point source, lower frequency, greater penetration
 - Trauma workhorse, cardio (fast response)
- Etc.....



Beam/inquiry patterns



Probe manipulation

- Probe marker orientation
 - Indicates location of specific beam
- Probe indicator positioning
 - Radiology convention: “Right hand generator rule”
 - Cardiology convention: “Left hand motor rule”
- Plane terminology:
 - Near field, far field
 - Leading edge, receding edge

You must think in 2D/plane terms!!

Probe manipulation


- This is important!
- **Fan** = tilt at same contact point
- **Sweep** = move into new plane
- **Slide** = move to extend same plane to side
- **Rock** = tilt to extend plane side to side
- **Rotate** = precise turn about central axis
 - PRACTICE!


Obtaining an image

- “Chance favors the prepared mind.” -Pasteur
- Mental preparation
 - What’s your question? What is relevant anatomy?
 - Select most appropriate probe
 - Apply gel (Benzoni’s rule: more gel)
 - Manipulate probe
 - Freeze or cine?

Artifacts

- Origins: timing, scatter, energy
- Timing creates **reverb**
 - Look for frequency/repetition effects
- Anisotropy
 - Beam **scatter**; look from several angles
- **Energy** absorption
 - Augmentation
 - Shadowing









Resources

- [Ohio State University – everything you need.](#)
- [U Virginia Medical U/S curriculum](#)
- [M-Turbo user manual](#)
- [Vascular access](#)
(Yes, nurses are doing this.)
- [Dawson/Mallen](#)
– Special ACEP offering; superb iBook.
- [Who can order a test? \(CMS\)](#)
- [AFP: POCUS in Family Medicine](#)
